Ventilátory ResMed a COVID-19

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Reference a poznámky

- ¹ Centers for Disease Control and Prevention. Interim Clinical Guidance for Management of Patients with Confirmed Coronavirus Disease (COVID-19). Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/hcp/clinical-guidance-management-patients.html> March 22, 2020
- ² WHO. Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected: interim guidance, 13 March 2020. Retrieved from <https://apps.who.int/iris/bitstream/handle/10665/331446/WHO-2019-nCoV-clinical-2020.4-eng.pdf?sequence=1&isAllowed=y> March 17, 2020.
- ³ Guan WJ, Ni ZY, Hu Y, et al. Clinical Characteristics of Coronavirus Disease 2019 in China [published online ahead of print, 2020 Feb 28]. *N Engl J Med.* 2020
- ⁴ Wang D, Hu B, Hu C, et al. Clinical Characteristics of 138 Hospitalized Patients With 2019 Novel Coronavirus-Infected Pneumonia in Wuhan, China [published online ahead of print, 2020 Feb 7]. *JAMA.* 2020
- ⁵ Murthy S, Gomersall CD, Fowler RA. Care for Critically Ill Patients With COVID-19 [published online ahead of print, 2020 Mar 11]. *JAMA.* 2020
- ⁶ NHS. Clinical guide for the use of acute non-invasive ventilation in adult patients hospitalised with suspected or confirmed coronavirus during the coronavirus pandemic [19 March 2020, Version 1].
- ⁷ Handbook of COVID-19 Prevention and Treatment. The First Affiliated Hospital, Zhejiang University School of Medicine Compiled According to Clinical Experience. Retrieved from <https://covid-19.alibabacloud.com/> on March 18, 2020.
- ⁸ Xie, J., Tong, Z., Guan, X. et al. Critical care crisis and some recommendations during the COVID-19 epidemic in China. *Intensive Care Med* (2020). <https://doi.org/10.1007/s00134-020-05979-7>
- ⁹ Medicines & Healthcare Products Regulatory Agency (MHRA). Rapidly Manufactured Ventilator Systems (RMVS). Retrieved from https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/874279/RMVS001_Rapidly_Manufactured_Ventilator_Specification__PDF.pdf March 23, 2020
- ¹⁰ US Department of Health and Human Services. Food and Drug Administration. Center for Devices and Radiological Health. Enforcement Policy for Ventilators and Accessories and Other Respiratory Devices During the Coronavirus Disease 2019 (COVID-19) Public Health Emergency. Guidance for Industry and Food and Drug Administration Staff. March 2020
- ¹¹ Hick, JL, Hanfling D, Wynia MK, et al. Duty to Plan: Health Care, Crisis Standards of Care, and Novel Coronavirus SARS-CoV-2. *NAM Perspectives*. Discussion paper. National Academy of Medicine. Washington, DC. 2020. Retrieved from <https://doi.org/10.31478/202003b> March 17, 2020.
- ¹² New York State Department of Health. Ventilator Allocation Guidelines. Retrieved from https://www.health.ny.gov/regulations/task_force/reports_publications/docs/ventilator_guidelines.pdf March 17, 2020.
- ¹³ Phillips Healthcare. Respiratory Management for Severe COVID-19 Patients Webinar – Part1. Retrieved from <https://www.youtube.com/watch?reload=9&v=rMT97toZNJw&feature=youtu.be> March 23, 2020.

¹⁴ Simonds AK, Hanak A, Chatwin M, et al. Evaluation of droplet dispersion during non-invasive ventilation, oxygen therapy, nebuliser treatment and chest physiotherapy in clinical practice: implications for management of pandemic influenza and other airborne infections. *Health Technol Assess.* 2010

¹⁵ Cascella M, Rajnik M, Cuomo A, et al. Features, Evaluation and Treatment Coronavirus (COVID-19) [Updated 2020 Mar 8]. In: *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing; 2020 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK554776/>

¹⁶ Gonzalez J, Maisonobe J, Oranger M, Mendoza-Ruiz A. Homecare respiratory equipment for patients suspected to be infected with the respiratory virus Covid19. Retrieved from <http://splf.fr/wp-content/uploads/2020/03/LES-PROCEDURES-DU-GAVO2-ProtectionVirale2020-MAJ14mars2020.pdf> March 18, 2020

¹⁷ Environmental Protection Agency. List N: disinfectants for use against SARS-CoV- 2. EPA website. Retrieved from <https://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2> March 22, 2020.

¹⁸ Australian Government Department of Health. Environmental cleaning and disinfecting principles – Version 1 (10/03/2020); Coronavirus disease (COVID-19). Retrieved from <https://www.health.gov.au/sites/default/files/documents/2020/03/environmental-cleaning-and-disinfection-principles-for-covid-19.pdf> Mar 18, 2020

¹⁹ Gonzalez J, Maisonobe J, Oranger M, Mendoza-Ruiz A. Homecare respiratory equipment for patients suspected to be infected with the respiratory virus Covid19. Retrieved from <http://splf.fr/wp-content/uploads/2020/03/LES-PROCEDURES-DU-GAVO2-ProtectionVirale2020-MAJ14mars2020.pdf> March 18, 2020

²⁰ CDC. Interim Infection Prevention and Control Recommendations for Patients with Suspected or Confirmed Coronavirus Disease 2019 (COVID-19) in Healthcare Settings. Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/infection-control/control-recommendations.html> March 17, 2020.

²¹ Environmental Protection Agency. List N: disinfectants for use against SARS-CoV-2. EPA website. Retrieved from <https://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2> March 22, 2020.

²² Australian Government Department of Health. Environmental cleaning and disinfecting principles – Version 1 (10/03/2020); Coronavirus disease (COVID-19).

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
ResMed

	Astral 100/150	Stellar 100/150	Lumis 150 ST-A	Lumis 150 ST
Patient Weight	>5 kg (all modes except iVAPS)	≥13 kg (all modes except iVAPS)	>13 kg (all modes except iVAPS)	>13 kg (all modes except iVAPS)
	>30 kg (iVAPS)	>30 kg (iVAPS)	>30 kg (iVAPS)	>30 kg (iVAPS)
Ventilation Delivery Specifications				
Patient Interface	non-invasive	non-invasive	non-invasive	non-invasive
	invasive	Invasive (with the use of ResMed leak valve)		
Circuit Type	Single valve	single leak	single leak	single leak
	Single leak			
	Double			
Pressure or Volume cycled	Pressure	Pressure	Pressure	Pressure
	Volume			
Operating Pressure Range	2-50 cmH2O: (S)T, P(A)C, iVAPS	2-40 cmH2O: S, ST, T, PAC, iVAPS	2 - 30 cmH2O: S, ST, T, PAC, iVAPS	2-25 cmH2O: S, ST, T, PAC, iVAPS
	3-20 cmH2O: CPAP	4-20 cmH2O: CPAP	4-20 cmH2O: CPAP	4-20 cmH2O: CPAP
	0-50 cmH2O: PS, PSIMV, V-SIMV			
Therapy Modes	CPAP, (S)T, P(A)C, (A)CV, P(A)CV, PSIMV, V-SIMV, PS, iVAPS	CPAP, S, ST (optional iBR),T, PAC, iVAPS	CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, ST (optional iBR),T, PAC, iVAPS
Performance Specifications				
Pressure Support Range		0-38 cmH2O	0-28 cmH2O	0-23 cmH2O
		0-30 cmH2O (in iVAPS mode)		
Intended Volume Range	Adult: 100 - 2,500 mL, Paed: 50 - 300 mL: (A)CV, V-SIMV	50-3000 mL	100–2500 mL	100–2500 mL
Supplementary Oxygen – Max Flow	30 L/min	30 L/min	15 L/min (S, ST, T, PAC, CPAP)	15 L/min (S, ST, T, PAC, CPAP)
		4 L/min (iVAPS)	4 L/min (iVAPS)	

	Astral 100/150	Stellar 100/150	Lumis 150 ST-A	Lumis 150 ST
Supplementary Oxygen – where to add	Device oxygen inlet	Device oxygen inlet	Device air outlet (eg. Side port of ClimateLineAir Oxy tube), mask	Device air outlet (eg. Side port of ClimateLineAir Oxy tube), mask
SpO2 measurement	Optional accessory	Optional accessory	Optional accessory	Optional accessory
FiO2 monitoring	Optional accessory	Optional accessory	n/a	n/a
Internal Battery Runtime	8 hours (a new battery under normal conditions)	8 hours (a new battery under normal conditions)	n/a	n/a
	Lithium-Ion battery, 14.4 V, 6.6 Ah, 95 Wh	Lithium-Ion battery, 14.4 V, 2.75 Ah, 40 Wh		
Humidification	External	Integrated or External	Integrated or External	Integrated or External
Safety				
Max single fault pressure	90 cmH2O	60 cm H2O (in all modes)	30 cmH2O for >6 s, or	30 cmH2O for >6 s, or
			40 cmH2O for >1 s	40 cmH2O for >1 s
Design Life	Device: 8 years	Device: 5 years	Device, PSU: 5 years	Device, PSU: 5 years
		Air Tubing: 6 months	Cleanable Humidifier: 2.5 years	Cleanable Humidifier: 2.5 years
			Air Tubing: 6 months	Air Tubing: 6 months
Alarms	Total power failure	Power fail	Power fail	n/a
	Circuit disconnection	Blocked tube	Blocked tube	
	Low Pressure, High Pressure	Tube disconnected, System fault	Tube disconnected	
	Obstruction	High temperature, Internal battery low/empty	High leak	

ResMed

	Astral 100/150	Stellar 100/150	Lumis 150 ST-A	Lumis 150 ST
	Apnea	Over pressure	Non-vented mask	
	Low / High MVe , Low / High MVi	High leak	Apnea	
	Low / High MVe , Low / High MVi	Non-vented mask	Low SpO2 (when oximeter connected)	
	Low / High Resp rate	Low minute ventilation	Low minute ventilation	
	High leak	Apnea		
	Ventilation stopped	High/Low pressure	System fault	
	Low / High SpO2	High/Low respiratory rate		
	Low / High pulse rate	High/Low FiO2		
	Low / High FiO2	Low SpO2		
	NV mask/Rebreathing			
	Incorrect circuit adapter			
	Critically low battery			
	Incorrect circuit attached			
	Safety reset complete			
	Battery inoperable			
	Low / High PEEP			
	Device overheating			
	Pressure line disconnected			
	Self-test failed			
	Flow sensor not calibrated			

	Astral 100/150	Stellar 100/150	Lumis 150 ST-A	Lumis 150 ST
	No SpO2 monitoring			
	No FiO2 monitoring			
	Internal battery degraded			
	Low internal battery			
	Circuit fault			
	Using internal battery			
	Battery 1 fault, Battery 2 fault			
	Power fault/No charging			
	PEEP blower failure			
	System Fault alarms			
Environmental				
Operating Temperature	32°F to 104°F (0°C to 40°C)	+32°F (0°C) to + 95°F (35°C)	+41°F to +95°F (+5°C to +35°C)	+41°F to +95°F (+5°C to +35°C)
Operating Humidity	5 to 93% noncondensing	10%–95% noncondensing	10 to 95% relative humidity, noncondensing	10 to 95% relative humidity, noncondensing
Storage and Transport Temperature	-13°F to 158°F (-25°C to 70°C)	-4°F (-20°C) to +140°F (60°C)	-4°F to +140°F (-20°C to +60°C)	-4°F to +140°F (-20°C to +60°C)
Storage and Transport Humidity	5 to 93% noncondensing	10%–95% noncondensing	5 to 95% relative humidity, noncondensing	5 to 95% relative humidity, noncondensing